Name -	
Partner	

# Viewing Cells under the Microscope

Purpose -	

### Procedure -

#### ONION CELL

- 1. Prepare a slide by placing two drops of iodine onto the slide.
- 2. Obtain a sample of onion. Snap the sample into <u>two</u> pieces and use tweezers to gently pull off a thin layer of onion skin. Place this sample onto a slide. Protect it with a cover slip.
- 3. Observe the onion sample and use a pencil draw your observations in the area below.

#### CHEEK CELL

- 1. Choose 1 group member to be the cell donor.
- 2. The donor obtains a sample by gently rubbing the inside of his/her cheek. Spread the sample onto a slide. Prepare the slide by placing <u>two</u> drops of iodine onto the area of the slide containing the sample. Protect your sample with a cover slip.
- 3. Observe the skin cells and use a pencil to draw your observations in the area below.

### Observations:

Onion Cell	Cell Type:
	Lens:
	Field of View:
	Size Calculation:

## <u>Cheek Cell</u>

	Cell Type:
	Lens:
	Field of View:
	Size Calculation:

### Analysis:

- 1. What parts of the cell were visible in your onion cell? Label these on your drawing.
- 2. What parts of the cell were visible in your cheek cell? Label these on your drawing.
- 3. What was the most obvious difference between the onion and cheek cells?
- 4. What type of cells are plant and animal cells classified as?

### Conclusion -